

SENDING SAMPLES TO GENOMEAST

Make sure to enclose your unique project ID (Sxxxx) in your shipment. Samples must not be sent before project validation by the Platform and by the Manager.

- **Samples must be registered in our interface** (<https://ngs-lims.igbmc.fr/>) **before sending**, including replacement samples.
- Please refer to our product sheets for volume and concentration requirements: <https://www.igbmc.fr/en/plateforms-and-services/platforms/genomeast#c3414>
- Please preferably label your tubes using the unique sample ID generated by our LIMS. Alternatively, you may use the same sample names registered in our interface.
- We recommend using 1.5 ml Eppendorf tubes for sample submission.
- For ready-to-load libraries as well as for low concentrated samples such as CHIP-seq samples, use « Low-bind » tubes.
- When sending frozen cell pellet, please schedule the shipment in advance with the platform so that the cells are frozen at -80 ° C as soon as they arrive on the platform.
- Dry ice is the recommended way to ship all samples and libraries except gDNA samples which must be sent with 4°C wet ice packs (blue ice).
- Ideally, next-day delivery services are preferred for most sample types. Favor **sending samples from Monday to Wednesday** so that the delivery date does not fall on a weekend.

All shipments should be sent to:

PLATEFORME GENOMEAST IGBMC
1 rue Laurent Fries
67404 Illkirch Cedex

In case of any queries, please contact us:

@ thibault@igbmc.fr
☎ (+33) 3 88 65 33 45